082207

English

RAPID SARS-COV-2 ANTIGEN TEST CARD
INSTRUCTION GUIDE FOR ANTERIOR NASAL SWAB SPECIMENS
For self-testing

	1N40C5-2	For 1 Test/Box
	1N40C5-4	For 5 Tests/Box
	1N40C5-6	For 20 Tests/Box

Please follow the instruction leaflet carefully.

INTENDED USE
Rapid SARS-CoV-2 Antigen Test Card is an immunochromatography based one step in vitro test. It is designed for the rapid qualitative determination of SARS-CoV-2 virus in anterior nasal swabs from individuals with or without symptoms or other epidemiological reasons to suspect COVID-19. Rapid SARS-CoV-2 Antigen Test Card detects the SARS-CoV-2 nucleocapsid protein (N protein). Rapid SARS-CoV-2 Antigen Test Card shall not be used as sole basis to diagnose or exclude SARS-CoV-2 infection. Children under 14 years of age should be assisted by an adult.

The novel coronavirus belongs to the β genus. COVID-19 is an acute respiratory infectious disease. People are generally susceptible. Currently, the patients infected by the novel coronavirus are the main source of infection, asymptomatic infected people can also be an infectious source. Based on the current epidemiological investigation, the incubation period is 1 to 14 days, mostly 3 to 7 days. The main manifestations include fever, fatigue and dry cough. Nasal congestion, runny nose, sore throat, myalgia and diarrhea are found in a few cases.

Components	For 1 Test/Box	For 5 Tests/Box	For 20 Tests/Box
Rapid SARS-CoV-2 Antigen Test Card (sealed foil pouch)	1	5	20
Sterile swab	1	5	20
Extraction tube	1	5	20
Sample extraction buffer	1	5	20
Tube stand (with its leaflet)	1	1	1
Tube stand (1 packaging)	1	1	1

PERFORMANCE (SENSITIVITY AND SPECIFICITY)
Rapid SARS-CoV-2 Antigen Test Card was compared to the confirmed clinical diagnosis. The study involved 1063 nasal samples. The test results are summarized below:

Evaluated Reagent Results	RT-PCR Results		Total
	Positive (+)	Negative (-)	
Positive (+)	425	1	426
Negative (-)	10	627	637
	435	628	1063

Sensitivity (PPA) = $425/435 \times 100\% = 97.70\%$ (95% CI: 96.29% - 99.11%)
Specificity (NPA) = $627/628 \times 100\% = 99.84\%$ (95% CI: 99.53% - 99.99%)
Accuracy (OPA) = $1052/1063 \times 100\% = 98.97\%$ (95% CI: 98.36% - 99.57%)

A feasibility study demonstrated that:
- 99.84% of non-professionals carried out the test without requiring assistance
- 99.82% of the different types of results were interpreted correctly

None of the following substances at the tested concentration showed any interference with the test.

Substance	Concentration
Whole Blood:	10%
Phenylephrine:	15%
Menthyl:	0.15%
Fluticasone Propionate:	5%
Oxetetracycline Phosphate:	0.5%
Alkaloid:	10%
Oxymetazolin:	0.0004%
Cromolyn:	5%
Mupirocin:	0.25%
Sodium chloride:	5%
Mucin:	2%
Oxymetazolin:	15%
Benzocaine:	0.15%
Zincam Nasal Spray:	5%
Human Anti-mouse Antibody (HAMA):	60 ng/mL

CROSS-REACTIVITY
Cross-reactivity of the Test Device was evaluated by testing viruses and other microorganisms. The final test concentrations of viruses and other microorganisms are documented in the Cross-Reactivity Study. The following viruses and other microorganisms except the Human SARS-coronavirus have no effect on the test results of the Test Device: Human coronavirus 229E, Human coronavirus OC43, Human coronavirus NL63, Human coronavirus HKU1, MERS coronavirus, Parainfluenza Virus 1-4, Enterovirus EV71, Respiratory syncytial virus, Rhinovirus, Influenza A virus (H1N1 and H3N2), Influenza B virus (Yamagata and Victoria), Adeno virus 71, Human Metapneumovirus (hMPV), Streptococcus pyogenes, Bordetella pertussis, Mycobacterium tuberculosis, Legionella pneumophila, Mycoplasma pneumoniae, Haemophilus influenzae, Candida albicans, Staphylococcus aureus, Pseudomonas aeruginosa, Escherichia coli, Pneumocystis jirovecii (PJP) and Pooled human nasal wash.

IMPORTANT INFORMATION BEFORE THE EXECUTION

1. Read this instruction guide carefully.
2. Do not use the product beyond the expiration date.
3. Do not use the test kit if the seal is damaged or the seal is broken.
4. Store the test device at 4 to 30 °C in the original sealed pouch. Do Not Freeze.
5. The product should be used at room temperature (15 °C to 30 °C). If the product has been stored in a cool area (less than 15 °C), leave it at normal room temperature for 30 minutes before use.
6. Handle all specimens as potentially infectious.
7. Inadequate or inappropriate specimen collection, storage, and transport may yield inaccurate test results.
8. Use the swabs included in the test kit to ensure optimal performance of the test.
9. Correct specimen collection is the most important step in the procedure. Make sure to collect enough specimen material (nasal secretion) with the swab, especially for anterior nasal sampling.
10. Do not use the swab several times before collecting specimen.
11. The specimens should be tested as soon as possible after collection.
12. Apply the drops of test specimen only to the specimen well (S).
13. Too many or too few drops of extraction solution can lead to an invalid or incorrect test result.
14. When used as intended, there should not be any contact with the extraction buffer. In case of contact with skin, eyes, mouth or other parts, rinse with clear water. If an irritation persists, consult a medical professional.
15. Children under 14 years of age should be assisted by an adult.

LIMITATIONS

1. The test is to be used exclusively for the qualitative detection of SARS-CoV-2 viral antigen in anterior nasal swab specimens. The exact concentration of SARS-CoV-2 viral antigen cannot be determined as part of this test.
2. Proper specimen collection is critical. Failure to follow the procedure may result in inaccurate test results. Improper collection, storage, or even freezing and thawing of the specimen can lead to inaccurate test results.
3. If the viral load of the specimen is below the detection limit of the test, the test may produce a negative result.
4. As with all diagnostic tests, a definitive clinical diagnosis should not be based on the result of a single test, but should be made by the physician after evaluation of all clinical and laboratory results.
5. A negative result does not exclude viral infection except for SARS-CoV-2 and should be confirmed by molecular diagnostic methods if COVID-19 is suspected.
6. A positive result does not exclude coinfection with other pathogens.
7. The SARS-CoV-2 rapid antigen test can detect both viable and non-viable SARS-CoV-2 material. The performance of the SARS-CoV-2 rapid test is dependent on viral load and other diagnostic methods performed on the same specimen.
8. Users should test specimens as soon as possible after specimen collection and within two hours of specimen collection.
9. Sensitivity for nasal or oropharyngeal swabs may be lower than nasopharyngeal swabs. It is recommended to use the nasopharyngeal swab specimens by healthcare professionals.
10. Monoclonal antibodies may fail to detect, or detect with less sensitivity, SARS-CoV-2 viruses that have undergone minor amino acid changes in the target epitope region.
11. The amount of antigen in a sample may decrease as the duration of illness increases. Viruses collected after day 5-7 of illness are more likely to be tested negative compared to a RT-PCR assay.
12. The test was validated with the assorted swabs. Use of alternative swabs may result in false negative results.
13. The validity of Rapid SARS-CoV-2 Antigen Test Card has not been proven for identification/confirmation of tissue culture isolates and should not be used in this capacity.
14. Positive test results do not rule out co-infections with other pathogens. Positive results may occur in cases of infection with SARS-CoV.

PREPARATION










- Clean, clean and dry a flat surface.
- Check the test kit contents. Make sure that nothing is damaged or broken.
- Rinse at hand.
- Blow your nose several times before collecting specimen.
- Wash hands.

DISPOSAL

The test kit may be disposed of with normal household waste in accordance with the applicable local regulations.

PROCEDURE

This test is suitable for people of all ages. The recommended ages are aging from 14 to 90. Children under 14 years of age should be tested by an adult. Do not continue the test if the child feels any pain.

**INTERPRETATION OF RESULTS**

Positive:
Two colored bands appear with one colored band in the Control Zone (C) and another in the Test Zone (T) within 15-20 minutes, the test result is positive. Caution: No matter how faint the colored band is in the Test Zone (T), the result should be considered as positive.

Negative:
If one colored band appears in the Control Zone (C) and no colored band appears in the Test Zone (T) within 15-20 minutes, the test result is negative.

Invalid:
If no color line appears in the control area (C) within 15-20 minutes, the test is invalid. Repeat the test with a new test card.

QUALITY CONTROL

The control line is an integrated reagent and is used to control the procedure. The control line appears when the test has been performed correctly and the reagents are reactive.

FREQUENTLY ASKED QUESTIONS (FAQ)

1. How does the detection work?
The reaction of the SARS-CoV-2 virus reacts with the stripe-like coating of the test line and, if present, results in a color change, i.e. a red line appears. Therefore, if the sample does not contain any viral proteins or antigens, there will be no red test line (T).

2. When should I test myself?
You can test yourself whether you have symptoms or not. Studies show that earlier testing within the first 4 days of illness typically means a higher viral load, which is easier to detect. Since the test result is a snapshot valid for that point in time, testing should be repeated as recommended by local authorities.

3. What can affect my test result? What should I pay attention to?
Be sure to blow your nose multiple times before collecting the specimen. Be sure to visibly collect sample material (nasal secretions). Perform the test immediately after taking the sample. Follow the instructions for use carefully. Apply the drops of extraction solution only to the sample well (S). Too many or too few drops of extraction solution can lead to an invalid or incorrect test result.

4. The test strip is clearly discolored or smudged? What is the reason for this?
Please note that the test card should not be used with more than 3 drops of sample, as the liquid absorption of the test strip is naturally limited. If the control line does not appear or the test strip is badly smudged or discolored, making it unreadable, please repeat the test according to the instructions.

5. I have taken the test, but I don't see a control line (C). What should I do?
Your test result is invalid. Observe the answer to question 4 and repeat the test according to the instructions for use.

6. I'm unsure about reading the result. What should I do?
For the result to be positive, 2 straight horizontal lines must be clearly visible with the full width of the cassette. If you are still unsure about the results, contact the nearest health authority according to the recommendations of your local authorities.

7. My result is positive. What should I do?
If your result is positive and the test kit clearly indicates the control line as well as the test line, you should contact the nearest medical facility as recommended by your local authorities or you should be double-checked and the authority or facility will explain the appropriate next steps.

8. My result is negative. What should I do?
If the test kit only clearly shows the control line, this may mean that you are negative or that the viral load is too low to be detected. If you experience symptoms (fever, fatigue, loss of sense of smell or taste, etc.), always consult your primary care physician, or the nearest health care facility as recommended by your local authorities. If you are not sure, you can repeat the test.

9. How can I dispose of the product?
The test kit may be disposed of with normal household waste in accordance with the applicable local regulations.

Accessories	Manufacturer	EC-Representative	CE-Mark
Swab A	Jiangsu Changfeng Medical Industry Co., Ltd. Touqiao Town, Guangming District Yangzhou 225109 Jiangsu P.R. China	Lins Service & Consulting GmbH Obere Seegasse 34/2,69124 Heidelberg Germany	 acc. 93/42/EEC
Swab B	Goodwood Medical Care Ltd. 1-2 Floor, 3-919 Yongzheng Street, Jinzhou District Dalian 116100 Liaoning China	CMC Medical Devices & Drugs S.L. C/ Horacio Lugo No18, CP 29006, Malaga, Spain	 acc. 93/42/EEC
Swab C	Zhejiang Gongdong Medical Technology Co., Ltd. No. 10 Beiyuan Ave., Huangyan 318020 Taizhou, Zhejiang, P.R.China	Shanghai International Holding Corp. GmbH (Europe) Eiffelstrasse 80, 20537 Hamburg, Germany	 acc. 93/42/EEC
Swab D	Jiangsu Hanheng Medical Technology Co., Ltd. 16-84#1 North Qingyang Road, Tianjing District, 213017, Changzhou, Jiangsu, China	Luxus Lebenswelt GmbH Kochstr. 1, 47877, Willich, Germany	 acc. 93/42/EEC

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ACCESORII:

Producător	Reprezentant european	Marcaj CE
Tampon A Jiangsu Changfeng Medical Industry Co., Ltd. Toujiao Town, Guangning District Yangzhou 225109 Jiangsu P.R. China	Lins Service & Consulting GmbH Obere Seegasse 34/2,69124 Heidelberg Germany	CE 1197 acc. 93/42/EEC
Tampon B Goodwood Medical Care Ltd. 1-2 Floor, 3-919 Yongzheng Street Jinzhou District Dalian 116100 Liaoning China	CMC Medical Devices & Drugs S.L. C/ Horacio Lengua No18, CP 29006, Málaga, Spain	CE 1197 acc. 93/42/EEC
Tampon C Zhejiang Gongdong Medical Technology Co., Ltd. No. 10 Beiuyan Ave., Huangyan 318020 Taizhou, Zhejiang, P.R. China	Shanghai International Holding Corp. GmbH (Europe) Eiffelstrasse 80, 20537 Hamburg, Germany	CE 1123 acc. 93/42/EEC
Tampon D Jiangsu Hanheng Medical Technology Co., Ltd. 16-B4-#1 North Qingyang Road, Tanning District, 213017, Changzhou, Jiangsu, China	Luxus Lebenswelt GmbH Kochstr. 1, 47877, Willich, Germany	CE 1197 acc. 93/42/EEC

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Producător: Xiamen Boson Biotech Co., Ltd.
Reprezentant autorizat: Lotus NL B.V.